



SLC38A1 mouse mAb

Catalog No	YP-mAb-17740
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC38A1
Alternative Names	ATA1; NAT2; SAT1; SNAT1
Research Field	Signal Transduction
Product Categories	Primary antibody
Host	mouse
Molecular Weight	Calculated MW: 54 kDa; Observed MW: 54 kDa
Clonality	Monoclonal Antibody
Clonality No.	-
Dilution	WB: 1/500-1/1000
Immunogen	Synthetic peptide of human SLC38A1
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	Functions as a sodium-dependent amino acid transporter.
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

