



GFAP mouse mAb

Catalog No	YP-mAb-17745
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GFAP
Immunogen	A synthesized peptide derived from human GFAP
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	glial fibrillary acidic protein;GFAP;FLJ45472;cb345;ALXDRD
Observed Band	Calculated MW: 50 kDa; Observed MW: 50 kDa
Calculated Molecular Weight	Calculated MW: 50 kDa; Observed MW: 50 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	GFAP is commonly used as a marker for intracranial and intraspinal tumors arising from astrocytes. In addition, GFAP intermediate filaments are also present in nonmyelin-forming Schwann cells in the peripheral nervous system
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

