



GIDRP88 Monoclonal Antibody

Catalog No	YP-mAb-01749
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	R3HCC1L
Immunogen	The antiserum was produced against synthesized peptide derived from human GIDRP88. AA range:721-770
Specificity	GIDRP88 Monoclonal Antibody detects endogenous levels of GIDRP88 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Putative mitochondrial space protein 32.1;PSORT;Growth inhibition and differentiation-related protein 88;GIDRP88;GIDRP86;R3HCC1L;R3H and coiled-coil domain-containing protein 1-like;C10orf28
Observed Band	88kD
Calculated Molecular Weight	88kD
Cell Pathway	exon-exon junction complex,
Tissue Specificity	Expressed in placenta.
Function	tissue specificity:Expressed in placenta.,
Background	tissue specificity:Expressed in placenta.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

