



Myf-5 Monoclonal Antibody

Catalog No	YP-mAb-01885
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MYF5
Immunogen	The antiserum was produced against synthesized peptide derived from human MYF5. AA range:21-70
Specificity	Myf-5 Monoclonal Antibody detects endogenous levels of Myf-5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Myf-5;MYF5;Myogenic factor 5;BHLHC2;Class C basic helix-loop-helix protein 2
Observed Band	28kD
Calculated Molecular Weight	28kD
Cell Pathway	Nucleus .
Tissue Specificity	Skeletal muscle,
Function	function:Involved in muscle differentiation (myogenic factor). Induces fibroblasts to differentiate into myoblasts. Probable sequence specific DNA-binding protein., similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein.,
Background	function:Involved in muscle differentiation (myogenic factor). Induces fibroblasts to differentiate into myoblasts. Probable sequence specific DNA-binding protein., similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein.,
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

