



ABHD8 Monoclonal Antibody

Catalog No	YP-mAb-02459
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ABHD8
Immunogen	The antiserum was produced against synthesized peptide derived from human ABHD8. AA range:291-340
Specificity	ABHD8 Monoclonal Antibody detects endogenous levels of ABHD8 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CGI58;IECN2;NCIE2;ABHD5;ABHD8;Abhydrolase domain-containing protein 8
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	extracellular exosome,
Tissue Specificity	Brain,Embryo,Lymph,Muscle,
Function	similarity:Belongs to the AB hydrolase superfamily.,
Background	This gene is upstream of, and in a head-to-head orientation with the gene for the mitochondrial ribosomal protein L34. The predicted protein contains alpha/beta hydrolase fold and secretory lipase domains. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

