



CYP7B1 Monoclonal Antibody

Catalog No	YP-mAb-02608
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CYP7B1
Immunogen	The antiserum was produced against synthesized peptide derived from human Cytochrome P450 7B1. AA range:101-150
Specificity	CYP7B1 Monoclonal Antibody detects endogenous levels of CYP7B1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	230160H5;24-hydroxycholesterol 7-alpha-hydroxylase;25/26-hydroxycholesterol 7-alpha-hydroxylase;3-hydroxysteroid 7-alpha hydroxylase;CBAS3;CYP7B1;25-hydroxycholesterol 7-alpha-hydroxylase;Cytochrome P450 7B1;Oxysterol 7-alpha-hydroxylase
Observed Band	58kD
Calculated Molecular Weight	58kD
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Microsome membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed. Expressed in brain, testis, ovary, prostate, liver, colon, kidney, small intestine, thymus and spleen.
Function	catalytic activity:Cholest-5-ene-3-beta,25-diol + NADPH + O(2) = cholest-5-ene-3-beta,7-alpha,25-triol + NADP(+) + H(2)O.,catalytic activity:Cholest-5-ene-3-beta,27-diol + NADPH + O(2) = cholest-5-ene-3-beta,7-alpha,27-triol + NADP(+) + H(2)O.,cofactor:Heme group.,disease:Defects in CYP7B1 are the cause of congenital bile acid synthesis defect type 3 (CBAS3) [MIM:603711]. Clinical features include severe cholestasis, cirrhosis and liver synthetic failure. Hepatic microsomal oxysterol 7-alpha-hydroxylase activity is undetectable., disease:Defects in CYP7B1 are the cause of spastic paraplegia autosomal



recessive type 5A (SPG5A) [MIM:270800]. Spastic paraplegia is a neurodegenerative disorder characterized by a slow, gradual, progressive weakness and spasticity of the lower limbs. Rate of progression and the severity of symptoms are quite variable. Initial symptoms may include difficulty with

Background

This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This endoplasmic reticulum membrane protein catalyzes the first reaction in the cholesterol catabolic pathway of extrahepatic tissues, which converts cholesterol to bile acids. This enzyme likely plays a minor role in total bile acid synthesis, but may also be involved in the development of atherosclerosis, neurosteroid metabolism and sex hormone synthesis. Mutations in this gene have been associated with hereditary spastic paraplegia (SPG5 or HSP), an autosomal recessive disorder. [provided by RefSeq, Apr 2016],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

