



GNPAT Monoclonal Antibody

Catalog No	YP-mAb-02646
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GNPAT
Immunogen	The antiserum was produced against synthesized peptide derived from human GNPAT. AA range:231-280
Specificity	GNPAT Monoclonal Antibody detects endogenous levels of GNPAT protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DAP AT;DAPAT;DHAP AT;DHAPAT;GNPAT;Dihydroxyacetone phosphate acyltransferase;DAP-AT;DHAP-AT;Acyl-CoA: dihydroxyacetonephosphateacyltransferase;Glycerone-phosphate O-acyltransferase
Observed Band	77kD
Calculated Molecular Weight	77kD
Cell Pathway	Peroxisome membrane ; Peripheral membrane protein ; Matrix side . Exclusively localized to the luminal side of the peroxisomal membrane. .
Tissue Specificity	Aorta endothelial cell,Brain,Liver,Lung,Thymus,
Function	catalytic activity:Acyl-CoA + glycerone phosphate = CoA + acylglycerone phosphate.,disease:Defects in GNPAT are the cause of rhizomelic chondrodysplasia punctata type 2 (RCDP2) [MIM:222765]. RDCP2 is characterized by rhizomelic shortening of femur and humerus, vertebral disorders, cataract, cutaneous lesions and severe mental retardation., domain:The HXXXXD motif is essential for acyltransferase activity and may constitute the binding site for the phosphate moiety of the glycerol-3-phosphate., pathway:Membrane lipid metabolism; glycerophospholipid metabolism., similarity:Belongs to the GPAT/DAPAT family.,subcellular location:Exclusively



localized to the luminal side of the peroxisomal membrane.,subunit:May be part of an heterotrimeric complex composed of DAP-AT, ADAP-S and a modified form of DAP-AT.,

Background

This gene encodes an enzyme located in the peroxisomal membrane which is essential to the synthesis of ether phospholipids. Mutations in this gene are associated with rhizomelic chondrodysplasia punctata. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2015],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

