



ARMCX3 Monoclonal Antibody

Catalog No	YP-mAb-03716
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ARMCX3
Immunogen	The antiserum was produced against synthesized peptide derived from human ARMX3. AA range:291-340
Specificity	ARMCX3 Monoclonal Antibody detects endogenous levels of ARMCX3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UNQ2517/PRO6007;Protein ALEX3;KIAA0443;BM-017;Armadillo repeat-containing X-linked protein 3;ARMCX3;ALEX3;ARM protein lost in epithelial cancers on chromosome X 3
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Mitochondrion outer membrane ; Single-pass membrane protein . Cytoplasm . Nucleus .
Tissue Specificity	Bone marrow,Brain,Epithelium,Kidney,Liver,Retina,Skin,Testis,
Function	similarity:Contains 3 ARM repeats.,
Background	armadillo repeat containing, X-linked 3(ARMCX3) Homo sapiens This gene encodes a member of the ALEX family of proteins which may play a role in tumor suppression. The encoded protein contains a potential N-terminal transmembrane domain and a single Armadillo (arm) repeat. Other proteins containing the arm repeat are involved in development, maintenance of tissue integrity, and tumorigenesis. This gene is closely localized with other family members on the X chromosome. Three transcript variants encoding the same



protein have been identified for this gene. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

