



DnaJC17 Monoclonal Antibody

Catalog No	YP-mAb-03817
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DNAJC17
Immunogen	The antiserum was produced against synthesized peptide derived from human DNAJC17. AA range:11-60
Specificity	DnaJC17 Monoclonal Antibody detects endogenous levels of DnaJC17 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DnaJ (Hsp40) homolog;subfamily C;member 17;DnaJC17;DnaJ homolog subfamily C member 17
Observed Band	35kD
Calculated Molecular Weight	35kD
Cell Pathway	Cytoplasm . Nucleus . Predominantly nuclear. .
Tissue Specificity	Brain,Epithelium,
Function	similarity:Contains 1 J domain.,similarity:Contains 1 RRM (RNA recognition motif) domain.,
Background	similarity:Contains 1 J domain.,similarity:Contains 1 RRM (RNA recognition motif) domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

