



GPR68 Monoclonal Antibody

Catalog No	YP-mAb-13347
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GPR68
Immunogen	The antiserum was produced against synthesized peptide derived from human SPR1. AA range:281-330
Specificity	GPR68 Monoclonal Antibody detects endogenous levels of GPR68 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	OGR-1;OGR1;OGR 1;G-protein coupled receptor 68;GPR12A;GPR68;Ovarian cancer G-protein coupled receptor 1;Sphingosylphosphorylcholine receptor
Observed Band	41kD
Calculated Molecular Weight	41kD
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	Found at low level in a wide range of tissues, but significantly expressed in lung, kidney, bone and nervous system.
Function	similarity:Belongs to the G-protein coupled receptor 1 family.,
Background	similarity:Belongs to the G-protein coupled receptor 1 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

