



VPAC1 Monoclonal Antibody

Catalog No	YP-mAb-13710
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	VIPR1
Immunogen	The antiserum was produced against synthesized peptide derived from human VIPR1. AA range:332-381
Specificity	VPAC1 Monoclonal Antibody detects endogenous levels of VPAC1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	II;HVR1;RDC1;V1RG;VIPR;VIRG;VAPC1;VPAC1;VPAC1R;VIP-R-1;VPCAP1R; PACAP-R2;PACAP-R-2;VIPR1;Vasoactive intestinal polypeptide receptor 1; Pituitary adenylate cyclase-activating polypeptide type II receptor;PACAP type II receptor
Observed Band	52kD
Calculated Molecular Weight	52kD
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	In lung, HT-29 colonic epithelial cells, Raji B-lymphoblasts. Lesser extent in brain, heart, kidney, liver and placenta. Not expressed in CD4+ or CD8+ T-cells. Expressed in the T-cell lines HARRIS, HuT 78, Jurkat and SUP-T1, but not in the T-cell lines Peer, MOLT-4, HSB and YT.
Function	function:This is a receptor for VIP. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase. The affinity is VIP = PACAP-27 > PACAP-38.,similarity:Belongs to the G-protein coupled receptor 2 family.,tissue specificity:In lung, HT29 colonic epithelial cells, Raji B-lymphoblasts. Lesser extent in brain, heart, kidney, liver and placenta. Not expressed in CD4+ or CD8+ T-cells. Expressed in the T-cell lines HARRIS, HuT 78, Jurkat and Tsup-1,



but not in the T-cell lines PEER, MOLT-4, HSB and YT.,

Background

vasoactive intestinal peptide receptor 1(VIPR1) Homo sapiens This gene encodes a receptor for vasoactive intestinal peptide, a small neuropeptide. Vasoactive intestinal peptide is involved in smooth muscle relaxation, exocrine and endocrine secretion, and water and ion flux in lung and intestinal epithelia. Its actions are effected through integral membrane receptors associated with a guanine nucleotide binding protein which activates adenylate cyclase. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

