



IRF2BP2 Mouse mAb

Catalog No	YP-mAb-19495
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	IRF2BP2
Immunogen	Synthesized peptide derived from human IRF2BP2. AA range:418-518
Specificity	This antibody detects endogenous levels of IRF2BP2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	IRF2BP2;Interferon regulatory factor 2-binding protein 2 (IRF-2-binding protein 2) (IRF-2BP2)
Observed Band	
Calculated Molecular Weight	65kD
Cell Pathway	Cytoplasm. Nucleus.
Tissue Specificity	
Function	Acts as a transcriptional corepressor in a IRF2-dependent manner; this repression is not mediated by histone deacetylase activities . Represses the NFAT1-dependent transactivation of NFAT-responsive promoters . Acts as a coactivator of VEGFA expression in cardiac and skeletal muscles . Plays a role in immature B-cell differentiation .
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images