



LMOD1 Mouse mAb

Catalog No	YP-mAb-19497
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	1D; D1; 64kD; MMIHS3; SMLMOD; SM-LMOD; LMOD1
Observed Band	
Calculated Molecular Weight	67kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	The leiomodlin 1 protein has a putative membrane-spanning region and 2 types of tandemly repeated blocks.The transcript is expressed in all tissues tested, with the highest levels in thyroid, eye muscle, skeletal muscle, and ovary. Increased expression of leiomodlin 1 may be linked to Graves' disease and thyroid-associated ophthalmopathy.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images