



VILIP3 Mouse mAb

Catalog No	YP-mAb-19502
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	HPCAL1 BDR1
Immunogen	Synthesized peptide derived from human VILIP3. AA range:1-56
Specificity	This antibody detects endogenous levels of VILIP3 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	HPCAL1 BDR1;Hippocalcin-like protein 1 (Calcium-binding protein BDR-1) (HLP2) (Visinin-like protein 3) (VILIP-3)
Observed Band	
Calculated Molecular Weight	21kD
Cell Pathway	Membrane ; Lipid-anchor .
Tissue Specificity	
Function	May be involved in the calcium-dependent regulation of rhodopsin phosphorylation.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

