




NCEH1 mouse mAb

Catalog No	YP-mAb-11477
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NCEH1 AADACL1 KIAA1363
Immunogen	Synthesized peptide derived from human NCEH1 AA range: 8-58
Specificity	This antibody detects endogenous levels of NCEH1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:3.1.1.71;EC:3.1.1.-;Arylacetamide deacetylase-like 1;Acetylalkylglycerol acetylhydrolase;AADACL 1;NCEH1;NCEH1 AADACL1 KIAA1363
Observed Band	
Calculated Molecular Weight	45kD
Cell Pathway	Cell membrane ; Single-pass type II membrane protein . Microsome .
Tissue Specificity	Expressed in monocyte-derived macrophages. Up-regulated in invasive melanoma and breast carcinoma cell lines.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using NCEH1 mouse mAb