



NADK Rabbit pAb

Catalog No	YP-Ab-18591
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	NADK
Protein Name	NAD kinase (Poly(P)/ATP NAD kinase)
Immunogen	Synthesized peptide derived from human NADK
Specificity	This antibody detects endogenous levels of NADK at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	49kD
Cell Pathway	
Tissue Specificity	Widely expressed but not detected in skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images