



Glutathione Reductase Mouse mAb

Catalog No	YP-mAb-17744
Isotype	IgG
Reactivity	Human, Mouse,Rat
Applications	WB
Gene Name	GSR
Alternative Names	GR;HEL-75;HEL-S-122m
Research Field	Signal Transduction
Product Categories	Primary antibody
Host	Rabbit
Molecular Weight	Calculated MW: 56 kDa; Observed MW: 51 kDa
Clonality	Monoclonal Antibody
Clonality No.	R09-5A9
Dilution	WB: 1/500-1/1000
Immunogen	A synthetic peptide of human Glutathione Reductase
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Form
Buffer System	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Concentration	1 mg/ml
Purity	≥90%
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	Maintains high levels of reduced glutathione in the cytosol.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

