

BMF Mouse mAb

货号	YP-mAb-18304
同位型	IgG
应用	WB
种属	Human,Mouse,Rat
基因名称	BMF
蛋白名称	Bcl-2-modifying factor
免疫原	Synthesized peptide derived from human BMF
特异性	This antibody detects endogenous levels of BMF at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Monoclonal,Mouse,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	20kD
功能	May play a role in apoptosis. Isoform 1 seems to be the main initiator.
组织表达	Isoform 1 is mainly expressed in B-lymphoid cells. Isoform 2 and isoform 3 are mainly expressed in B-CLL and normal B-cells.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images