



Collagen Type I A1(ABT162) Mouse mAb

货号	YP-Ab-17738
同位型	IgG
应用	IHC;WB;IF
种属	Human;Mouse;Rat
靶点	Collagen I
简介	>>PI3K-Akt signaling pathway;>>Focal adhesion;>>ECM-receptor interaction; >>Platelet activation;>>Relaxin signaling pathway;>>AGE-RAGE signaling pathway in diabetic complications;>>Protein digestion and absorption; >>Amoebiasis;>>Human papillomavirus infection;>>Proteoglycans in cancer; >>Diabetic cardiomyopathy
基因名称	COL1A1
免疫原	Synthesized peptide derived from human Collagen Type I AA range: 1200-1464
特异性	The antibody can specifically recognize human Collagen Type I protein, collagen types II, III, IV and V do not respond to the antibody.
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Mouse, Monoclonal/IgG2b, Kappa
稀释	IHC-p 1:100-500, WB 1:200-1000, IF 1:100-500
纯化工艺	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
其他名称	COL1A1;1E9A7;Alpha 1 type I collagen;Alpha-1 type I collagen;COL I;Collagen I; Collagen Type I;Collagen alpha-1(I) chain (Alpha-1 type I collagen)
分子量	139 kDa
背景	This gene encodes the pro-alpha1 chains of type I collagen whose triple helix comprises two alpha1 chains and one alpha2 chain. Type I is a fibril-forming collagen found in most connective tissues and is abundant in bone, cornea, dermis and tendon. Mutations in this gene are associated with osteogenesis imperfecta types I-IV, Ehlers-Danlos syndrome type VIIA, Ehlers-Danlos syndrome Classical type, Caffey Disease and idiopathic osteoporosis. Reciprocal translocations between chromosomes 17 and 22, where this gene and the gene for platelet-derived growth factor beta are located, are associated with a particular type of skin tumor called dermatofibrosarcoma protuberans, resulting from unregulated expression of the growth factor. Two transcripts, resulting from the use of alternate polyadenylation signals, have been identified for this gene. [provided by R. Dalgleish, Feb 2008],
功能	disease:A chromosomal aberration involving COL1A1 is a cause of dermatofibrosarcoma protuberans (DFSP) [MIM:607907]. Translocation t(17;22)(q22;q13) with PDGF. DFSP is an uncommon, locally aggressive, but rarely metastasizing tumor of the deep dermis and subcutaneous tissue. It typically



occurs during early or middle adult life and is most frequently located on the trunk and proximal extremities.,disease:Defects in COL1A1 are a cause of Ehlers-Danlos syndrome type 1 (EDS1) [MIM:130000]; also known as Ehlers-Danlos syndrome gravis. EDS is a connective tissue disorder characterized by hyperextensible skin, atrophic cutaneous scars due to tissue fragility and joint hyperlaxity. EDS1 is the severe form of classic Ehlers-Danlos syndrome., disease:Defects in COL1A1 are a cause of osteogenesis imperfecta type I (OI-I) [MIM:166200]. OI-I is a dominantly inherited serious newborn disease character

细胞定位	Secreted, extracellular space, extracellular matrix .
组织表达	Forms the fibrils of tendon, ligaments and bones. In bones the fibrils are mineralized with calcium hydroxyapatite.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



