



TRMT10A rabbit pAb

货号	YP-Ab-17739
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
靶点	TRMT10A
基因名称	TRMT10A RG9MTD2
免疫原	Synthesized peptide derived from human TRMT10A
特异性	This antibody detects endogenous levels of TRMT10A at Human, Mouse,Rat
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	339 aa, 40 kDa
功能	S-adenosyl-L-methionine-dependent guanine N(1)-methyltransferase that catalyzes the formation of N(1)-methylguanine at position 9 (m1G9) in tRNAs . Probably not able to catalyze formation of N(1)-methyladenine at position 9 (m1A9) in tRNAs .
细胞定位	Nucleus . Nucleus, nucleolus .
组织表达	Expressed in embryonic and fetal brain. It is expressed throughout the dorsal telencephalon at 8 and 11 weeks of gestation, with highest expression in ventricular zone and marginal zone. Detected in cerebellar cortex and nuclei, but not in dorsal telencephalon, at later stages.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

