



PCYT2 mouse mAb

Catalog No	YP-mAb-17123	
Isotype	IgG	
Reactivity	Human;Mouse;Rat	
Applications	WB	
Gene Name	type	Monoclonal Antibody
Immunogen	Recombinant fusion protein of human PCYT2 (NP_002852.1).	
Specificity	matters needing attention	Avoid repeated freezing and thawing!
Formulation		
Source	Host	Mouse
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.	
Dilution	WB 1:500-1:2000	
Concentration	1 mg/ml	
Purity	≥90%	
Storage Stability	-20°C/1 year	
Synonyms	type	Monoclonal Antibody
Observed Band	Detects a band of approximately 44kD	
Calculated Molecular Weight	Detects a band of approximately 44kD	
Cell Pathway		
Tissue Specificity	Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.
Function	Disclaimer note	The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors.
Background		
matters needing attention	Avoid repeated freezing and thawing!	

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images