



mouse anti-CEPT1 (N-term) Monoclonal Antibody

Catalog No	YP-mAb-17124	
Isotype	IgG	
Reactivity	Human;Mouse;Rat	
Applications	WB	
Gene Name	Predict the reaction species	Mouse, Rat, Xenopus
Immunogen		
Specificity	This CEPT1 antibody is generated from Mouses immunized with a KLH conjugated synthetic peptide between 29-57 amino acids from the N-terminal region of human CEPT1.	
Formulation		
Source	Host	Mouse
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.	
Dilution	WB 1:500-1:2000	
Concentration	1 mg/ml	
Purity	≥90%	
Storage Stability	-20°C/1 year	
Synonyms	Predict the reaction species	Mouse; Rat; Xenopus
Observed Band	Predicted molecular weight: 47kD	
Calculated Molecular Weight	19kD	
Cell Pathway	Endoplasmic reticulum membrane; Multi-pass membrane protein. Nucleus membrane; Multi-pass membrane protein	
Tissue Specificity	matters needing attention	Avoid repeated freezing and thawing!
Function	Disclaimer note The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors.	
Background		
matters needing attention	Avoid repeated freezing and thawing!	

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images