



VEGF Receptor 1 mouse mAb

Catalog No	YP-mAb-17140	
Isotype	IgG	
Reactivity	Human;Mouse;Rat	
Applications	WB	
Gene Name	FLT1	
Immunogen	A synthesized peptide derived from human VEGF Receptor 1	
Specificity		
Formulation		
Source		
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.	
Dilution	WB 1:500-1:2000	
Concentration	1 mg/ml	
Purity	≥90%	
Storage Stability	-20°C/1 year	
Synonyms	FLT1;Alternative Names VEGFR-1;VEGF Receptor 1;FLT-1;FLT;FRT; VEGFR1;Fms-like tyrosine kinase 1;Tyrosine-protein kinase FRT; Modification Unmodified	
Observed Band	Form	Liquid
Calculated Molecular Weight	Form	Liquid
Cell Pathway		
Tissue Specificity	Buffer System Mouse IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.	
Function		
Background	Gene ID:2321	
matters needing attention	Avoid repeated freezing and thawing!	
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.	



Products Images