



Gz-α (Phospho-Ser16) rabbit pAb

Catalog No	YP-Ab-10536
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	GNAZ
Immunogen	Synthesized peptide derived from human Gz-α (Phospho-Ser16)
Specificity	This antibody detects endogenous levels of Gz-α (Phospho-Ser16) at Human, Mouse, Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.192% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Gz-α; GNAZ; Gz-α (Phospho-Ser16); Guanine nucleotide-binding protein G(z) subunit alpha (G(x) alpha chain) (Gz-alpha)
Observed Band	41kD
Calculated Molecular Weight	41kD
Cell Pathway	Membrane; Lipid-anchor.
Tissue Specificity	Brain, PNS,
Function	function: Guanine nucleotide-binding proteins (G proteins) are involved as modulators or transducers in various transmembrane signaling systems., similarity: Belongs to the G-alpha family. G(i/o/t/z) subfamily., subunit: G proteins are composed of 3 units; alpha, beta and gamma. The alpha chain contains the guanine nucleotide binding site.,
Background	G protein subunit alpha z (GNAZ) Homo sapiens The protein encoded by this gene is a member of a G protein subfamily that mediates signal transduction in pertussis toxin-insensitive systems. This encoded protein may play a role in maintaining the ionic balance of perilymphatic and endolymphatic cochlear fluids. [provided by RefSeq, Jul 2008],


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of SY-SY5Y using primary antibody at 1:1000 dilution 4°C, overnight. Secondary antibody(catalog#:RS23920) was diluted at 1:10000 25 °C, 1.5hours