




ISK7 rabbit pAb

Catalog No	YP-Ab-12564
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	SPINK7 ECG2 UNQ745/PRO1474
Immunogen	Synthesized peptide derived from human ISK7 AA range: 1-50
Specificity	This antibody detects endogenous levels of ISK7 at Human, Mouse, Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.55% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ECG2; ECRG 2; ECRG2; ECRG-2; SPINK7; ISK7; SPINK7 ECG2 UNQ745/PRO1474; Serine protease inhibitor Kazal-type 7 (Esophagus cancer-related gene 2 protein) (ECRG-2)
Observed Band	20kD
Calculated Molecular Weight	20kD
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

hela 293T



Western Blot analysis of hela 293T using primary antibody at 1:1000 dilution 4°C, overnight. Secondary antibody(catalog#:RS23920) was diluted at 1:10000 25 °C, 1.5hours