



eIF2α (phospho Ser51) Polyclonal Antibody

Catalog No	YP-Ab-03515
Isotype	IgG
Reactivity	Human;Mouse;Rat;Dog;Pig;Fish
Applications	IF;WB;IHC;ELISA
Gene Name	EIF2S1
Immunogen	The antiserum was produced against synthesized peptide derived from human eIF2 alpha around the phosphorylation site of Ser51. AA range:21-70
Specificity	Phospho-eIF2α (S51) Polyclonal Antibody detects endogenous levels of eIF2α protein only when phosphorylated at S51.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IF: 1:50-200 Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EIF2S1/EIF2A;EIF 2;eIF 2 alpha;EIF 2A;EIF 2alpha;eIF2α;EIF2S1;Eukaryotic translation initiation factor 2 subunit 1;EIF2A;Eukaryotic translation initiation factor 2 subunit alpha;eIF-2-alpha;eIF-2A;eIF-2alpha
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Cytoplasm, Stress granule . Colocalizes with NANOS3 in the stress granules. .
Tissue Specificity	B cells,Brain,Fibroblast,Placenta,
Function	function:Functions in the early steps of protein synthesis by forming a ternary complex with GTP and initiator tRNA. This complex binds to a 40S ribosomal subunit, followed by mRNA binding to form a 43S preinitiation complex. Junction of the 60S ribosomal subunit to form the 80S initiation complex is preceded by hydrolysis of the GTP bound to eIF-2 and release of an eIF-2-GDP binary complex. In order for eIF-2 to recycle and catalyze another round of initiation, the GDP bound to eIF-2 must exchange with GTP by way of a reaction catalyzed by eIF-2B.,PTM:Substrate for at least 4 kinases: EIF2AK3/PERK, GCN2, HRI and PKR. Phosphorylation stabilizes the eIF-2/GDP/eIF-2B complex



and prevents GDP/GTP exchange reaction, thus impairing the recycling of eIF-2 between successive rounds of initiation and leading to global inhibition of translation. In case of infection by vaccinia virus or rotavirus

Background

The translation initiation factor EIF2 catalyzes the first regulated step of protein synthesis initiation, promoting the binding of the initiator tRNA to 40S ribosomal subunits. Binding occurs as a ternary complex of methionyl-tRNA, EIF2, and GTP. EIF2 is composed of 3 nonidentical subunits, the 36-kD EIF2-alpha subunit (EIF2S1), the 38-kD EIF2-beta subunit (EIF2S2; MIM 603908), and the 52-kD EIF2-gamma subunit (EIF2S3; MIM 300161). The rate of formation of the ternary complex is modulated by the phosphorylation state of EIF2-alpha (Ernst et al., 1987 [PubMed 2948954]).[supplied by OMIM, Feb 2010],

matters needing attention


Avoid repeated freezing and thawing!

Usage suggestions


This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images




Immunofluorescence analysis of A549. 1, primary Antibody (red) was diluted at 1:200 (4°C overnight). 2, Goat Anti Rabbit IgG (H&L) - Alexa Fluor 594 Secondary antibody was diluted at 1:1000 (room temperature, 50min). 3, Picture B: DAPI (blue) 10min.




Lei, Ting, et al. "Stress kinases, endoplasmic reticulum stress, and Alzheimer's disease related markers in peripheral blood mononuclear cells from subjects with increased body weight." *Scientific reports* 6 (2016): 30890.



Chen, Yong, et al. "Overexpression of long non coding RNA NORAD promotes invasion and migration in malignant melanoma via regulating the MIR 205 - ERLN2 pathway." *Cancer medicine* (2019).



Zheng, Qing, et al. "Isorhamnetin protects against bleomycin-induced pulmonary fibrosis by inhibiting endoplasmic reticulum stress and epithelial-mesenchymal transition." *International journal of molecular medicine* 43.1 (2019): 117-126.



Western blot analysis of lysates from 1) KB, 2) MCF-7 cells, (Green) primary antibody was diluted at 1:1000, 4 $^{\circ}$ over night, Dylight 800 secondary antibody(Immunoway:RS23920)was diluted at 1:10000, 37 $^{\circ}$ 1hour. (Red) Tubulin β Monoclonal Antibody(5G3) (Immunoway:YM3030) antibody was diluted at 1:5000 as loading control, 4 $^{\circ}$ over night,Dylight 680 secondary antibody(Immunoway:RS23710)was diluted at 1:10000, 37 $^{\circ}$ 1hour.