




PADI4 Mouse mAb

货号	YP-mAb-18313
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
免疫原	Recombinant fusion protein containing a sequence corresponding to amino acids 290-663 of human PADI4
来源	Monoclonal,Mouse,IgG
稀释	WB 1:500 - 1:2000
纯化工艺	Affinity purification
分子量	663 aa, 74 kDa
背景	This gene is a member of a gene family which encodes enzymes responsible for the conversion of arginine residues to citrulline residues. This gene may play a role in granulocyte and macrophage development leading to inflammation and immune response.
浓度	1 mg/ml
储存	StorageStability-20°C/1year
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of various lysates using PADI4 Mouse pAb (A16188) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane.

Blocking bufer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90s.