




RACK1 mouse mAb

货号	YP-mAb-18321
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
其他名称	Gnb2 rs1;GNB2L1;H12.3;HLC 7;PIG21;RACK1;HLC-7;Gnb2-rs1
免疫原	Recombinant fusion protein containing a sequence corresponding to amino acids 220-310 of human RACK1 (P63244)
稀释	WB 1:1000 - 1:6000
纯化工艺	Affinity purification
分子量	36 kDa
背景	Enables several functions, including cyclin binding activity; enzyme binding activity; and protein domain specific binding activity. Involved in several processes, including positive regulation of hydrolase activity; regulation of cellular protein metabolic process; and regulation of signal transduction. Located in several cellular components, including midbody; perinuclear region of cytoplasm; and phagocytic cup. Part of IRE1-RACK1-PP2A complex
细胞定位	
组织表达	
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of various lysates using RACK1 Mouse mAb (A3808) at 1:1000 dilution.


Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25 μ g per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 1s.



Western blot analysis of various lysates using RACK1 Mouse mAb (A3808) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25 μ g per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 10s.