


HECTD1 mouse mAb

货号	YP-mAb-18328
同位型	IgG
应用	WB
种属	Human;Mouse;Rat
其他名称	E3 ligase for inhibin receptor;E3 ubiquitin-protein ligase HECTD1;EC:2.3.2.26; EULIR;HECT domain containing 1;HECTD1
免疫原	Recombinant fusion protein containing a sequence corresponding to amino acids 2518-2610 of human HECTD1 (NP_056197.3)
稀释	WB 1:500-2000
纯化工艺	Affinity purification
分子量	2612 aa, 290 kDa
背景	Predicted to enable ubiquitin protein ligase activity. Predicted to be involved in anatomical structure development; positive regulation of proteasomal ubiquitin-dependent protein catabolic process; and protein K63-linked ubiquitination. Predicted to act upstream of or within several processes, including animal organ development; negative regulation of protein localization to plasma membrane; and protein autoubiquitination.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of various lysates using HECTD1 Mouse pAb (A9433) at 1:3000 dilution.

Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25μg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Enhanced Kit (RM00021).

Exposure time: 90s.