



MARCH1 Monoclonal Antibody

Catalog No	YP-mAb-02668
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	1-3
Immunogen	Synthesized peptide derived from the Internal region of human 40603.
Specificity	MARCH1 Monoclonal Antibody detects endogenous levels of 40603 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	44986;E3 ubiquitin-protein ligase MARCH1;MARCH1;RNF171;Membrane-associated RING finger protein 1;Membrane-associated RING-CH protein I; MARCH-I;RING finger protein 171
Observed Band	32kD
Calculated Molecular Weight	32kD
Cell Pathway	Golgi apparatus, trans-Golgi network membrane ; Multi-pass membrane protein . Lysosome membrane ; Multi-pass membrane protein . Cytoplasmic vesicle membrane ; Multi-pass membrane protein . Late endosome membrane ; Multi-pass membrane protein . Early endosome membrane ; Multi-pass membrane protein . Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in antigen presenting cells, APCs, located in lymph nodes and spleen. Also expressed in lung. Expression is high in follicular B-cells, moderate in dendritic cells and low in splenic T-cells.
Function	domain:The RING-CH-type zinc finger domain is required for E3 ligase activity., function:E3 ubiquitin-protein ligase that may mediate ubiquitination of TFRC, CD86 and FAS, and promote their subsequent endocytosis and sorting to lysosomes via multivesicular bodies. E3 ubiquitin ligases accept ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly



transfer the ubiquitin to targeted substrates.,pathway:Protein modification;
protein ubiquitination.,similarity:Contains 1 RING-CH-type zinc finger.,tissue
specificity:Expressed in lymph nodes, spleen and lung.,

Background

MARCH1 is a member of the MARCH family of membrane-bound E3 ubiquitin ligases (EC 6.3.2.19). MARCH proteins add ubiquitin (see MIM 191339) to target lysines in substrate proteins, thereby signaling their vesicular transport between membrane compartments. MARCH1 downregulates the surface expression of major histocompatibility complex (MHC) class II molecules (see MIM 142880) and other glycoproteins by directing them to the late endosomal/lysosomal compartment (Bartee et al., 2004 [PubMed 14722266]; Thibodeau et al., 2008 [PubMed 18389477]; De Gassart et al., 2008 [PubMed 18305173]) .[supplied by OMIM, Mar 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using MARCH1 Monoclonal Antibody