



UFL1 Rabbit pAb

Catalog No	YP-Ab-19196
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	UFL1 KIAA0776 NLBP
Immunogen	Synthesized peptide derived from human UFL1
Specificity	This antibody detects endogenous levels of UFL1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NLBP;RCAD;Maxer;KIAA0776;UFL1;UFL1 KIAA0776 NLBP;E3 UFM1-protein ligase 1 (LZAP-binding protein)
Observed Band	90kDa
Calculated Molecular Weight	90kDa
Cell Pathway	Endoplasmic reticulum membrane . Cytoplasm, cytosol . Nucleus . Chromosome . Recruited to double-strand breaks by the MRE11-RAD50-NBN (MRN) complex following DNA damage. .
Tissue Specificity	Ubiquitously expressed, with a high expression in liver (at protein level) (PubMed:20018847). Low expression in several invasive hepatocellular carcinomas, such Hep-G2, Hep 3B2.1-7, HLE and PLC (PubMed:20018847).
Function	E3 protein ligase that mediates ufmylation, the covalent attachment of the ubiquitin-like modifier UFM1 to lysine residues on target proteins, and which plays a key role in reticulophagy (also called ER-phagy) induced in response to endoplasmic reticulum stress . In response to endoplasmic reticulum stress, recruited to the endoplasmic reticulum membrane by DDRGK1, and mediates ufmylation of proteins such as RPN1 and RPL26/uL24, thereby promoting reticulophagy of endoplasmic reticulum sheets . Ufmylation-dependent reticulophagy inhibits the unfolded protein response (UPR) via ERN1/IRE1-



alpha . Ufmylation in response to endoplasmic reticulum stress is essential for processes such as hematopoiesis, blood vessel morphogenesis or inflammatory response . Regulates inflammation in response to endoplasmic reticulum stress by promoting reticulophagy, leading to inhibit the activity of the NF-kappa-B transcription factor (By similarity). Mediates ufmylation of DDRGK1 and CDK5RAP3; the role of these modifications is however unclear: as both DDRGK1 and CDK5RAP3 act as substrate adapters for ufmylation, it is uncertain whether ufmylation of these proteins is a collateral effect or is required for ufmylation . Catalyzes ufmylation of various subunits of the ribosomal complex or associated components, such as RPS3/uS3, RPS20/uS10, RPL10/uL16, RPL26/uL24 and EIF6 (By similarity). Anchors CDK5RAP3 in the cytoplasm, preventing its translocation to the nucleus which allows expression of the CCND1 cyclin and progression of cells through the G1/S transition . Also involved in the response to DNA damage: recruited to double-strand break sites following DNA damage and mediates monoufmylation of histone H4 . Catalyzes ufmylation of TRIP4, thereby playing a role in nuclear receptor-mediated transcription . Required for hematopoietic stem cell function and hematopoiesis (By similarity). Required for cardiac homeostasis (By similarity).

Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

