



Chimaerin 2 Rabbit pAb

Catalog No	YP-Ab-19507
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	CHN2 ARHGAP3 BCH
Immunogen	Synthesized peptide derived from human Chimaerin 2. AA range:1-130
Specificity	This antibody detects endogenous levels of Chimaerin 2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	CHN2 ARHGAP3 BCH;Beta-chimaerin (Beta-chimerin) (Rho GTPase-activating protein 3)
Observed Band	
Calculated Molecular Weight	51kD
Cell Pathway	Membrane ; Peripheral membrane protein .
Tissue Specificity	Highest levels in the brain and pancreas. Also expressed in the heart , placenta , and weakly in the kidney and liver. Expression is much reduced in the malignant gliomas , compared to normal brain or low-grade astrocytomas.
Function	GTPase-activating protein for p21-rac. Insufficient expression of beta-2 chimaerin is expected to lead to higher Rac activity and could therefore play a role in the progression from low-grade to high-grade tumors.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images