



TFF3 Mouse mAb

Catalog No	YP-mAb-19508
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	ITF; P1B;TFI;TFF3
Observed Band	
Calculated Molecular Weight	9kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Members of the trefoil family are characterized by having at least one copy of the trefoil motif, a 40-amino acid domain that contains three conserved disulfide bonds.They are stable secretory proteins expressed in gastrointestinal mucosa.Their functions are not defined, but they may protect the mucosa from insults, stabilize the mucus layer and affect healing of the epithelium.This gene is expressed in goblet cells of the intestines and colon.This gene and two other related trefoil family member genes are found in a cluster on chromosome 21.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images