



NPM3 Mouse mAb

Catalog No	YP-mAb-19509
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	NPM3
Immunogen	Synthesized peptide derived from human NPM3. AA range:40-110
Specificity	This antibody detects endogenous levels of NPM3 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	NPM3;Nucleoplasmin-3
Observed Band	
Calculated Molecular Weight	20kD
Cell Pathway	Nucleus . Nucleus , nucleolus . Mainly found in the granular component of the nucleolus. .
Tissue Specificity	Ubiquitous.
Function	Plays a role in the regulation of diverse cellular processes such as ribosome biogenesis , chromatin remodeling or protein chaperoning . Modulates the histone chaperone function and the RNA-binding activity of nucleolar phosphoprotein B23/NPM . Efficiently mediates chromatin remodeling when included in a pentamer containing NPM3 and NPM .
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images