



PHC2 mouse mAb

Catalog No	YP-mAb-10158
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PHC2 EDR2 PH2
Immunogen	Synthesized peptide derived from human PHC2 AA range: 233-283
Specificity	This antibody detects endogenous levels of human PHC2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EDR2;HPH2;PH2;Polyhomeotic like protein 2;Polyhomeotic-like protein 2;PHC2; PHC2 EDR2 PH2;Polyhomeotic-like protein 2 (hPH2) (Early development regulatory protein 2)
Observed Band	40 kDa
Calculated Molecular Weight	858 aa, 91 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

