



优品生物
UpingBio

网址: www.upingBio.com

邮箱: Upingbio.163.com

技术电话: 400-999-8863

小鼠抗猪免疫球蛋白G (IgG) 抗体 (生物素标记) 说明书

首先感谢您选择优品生物, 我司为您提供优质的免疫学试剂, 在您 使用本产品前, 请认真阅读说明书, 不清楚之处, 请联系技术人员解答, 这对您顺利完成实验很有帮助。

名 称: 小鼠抗猪免疫球蛋白G (IgG) 抗体 (生物素标记)

目 录 号: YPE800095-S

亚 型: IgG

靶 点: Pig IgG

组 成: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

来 源: Monoclonal, Mouse

稀 释 : WB 1:1000-10000, ELISA 1:10000-20000, IHC 1:50-100 IF 1:50-1000

纯化工艺: The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.

其他名称: BIOTIN*Monoclonal mouse Anti-Pig IgG

背 景 : Immunoway secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application. Serum adsorbed secondary antibodies are also available and are recommended for use with immunoglobulin-rich samples.

应用范围: WB;IHC;IF;ELISA

种 属: Pig

产品纯度: 2mg/mL。

储存方法: -15°C to -25°C/1 year(Do not lower than -25°C)

注意事项: Avoid repeated freezing and melting。

有效期限: Please use it within the validity period of the pipe body。

使用建议: This product can be used for immunological response related experiments, for more information, please consult technicians。