



优品生物
UpingBio

网址: www.upingBio.com

邮箱: Upingbio.163.com

技术电话: 400-999-8863

山羊抗小鼠免疫球蛋白 G (IgG) 抗体 (H+L) (生物素标记)

说明书

首先感谢您选择优品生物，我司为您提供优质的免疫学试剂，在您使用本产品前，请认真阅读说明书，不清楚之处，请联系技术人员解答，这对您顺利完成实验很有帮助。

名 称：山羊抗小鼠免疫球蛋白 G (IgG) 抗体 (H+L) (生物素标记)。

目 录 号：YP848536-S

亚 型：IgG

应 用：WB;IHC;IF;ELISA

种 属：Mouse

靶 点：Mouse IgG(H+L)

组 成：Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

来 源：Goat

稀 释：WB 1:1000-10000, ELISA 1:10000-20000, IHC 1:50-100 IF 1:50-1000

纯化工艺：The antibody was affinity-purified from goat antiserum by affinity-chromatography using epitope-specific immunogen.

浓 度：4 mg/ml

储 存：-15° C to -25° C/1 year(Do not lower than -25° C)

其他名称：BIOTIN* Goat Anti-Mouse IgG(H+L)

背 景：Immunoway secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application. Serum adsorbed secondary antibodies are also available and are recommended for use with immunoglobulin-rich samples.