



COAA1 Polyclonal Antibody

Catalog No	YP-Ab-07735
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	COL10A1
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	COAA1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Col10;Col10a-1;COAA1_BOVIN;COL10A1;A1L4P2_HUMAN;COAA1_HUMAN;collagen type X alpha 1 chain;collagen;type X;alpha 1;Schmid metaphyseal chondrodysplasia;COAA1;Collagen alpha-1(X) chain
Observed Band	66kDa
Calculated Molecular Weight	66kDa
Cell Pathway	Secreted, extracellular space, extracellular matrix .
Tissue Specificity	
Function	disease:Defects in COL10A1 are the cause of Schmid type metaphyseal chondrodysplasia (SMCD) [MIM:156500]. SMCD is a dominantly inherited disorder of the osseous skeleton. The cardinal features of the phenotype are mild short stature, coxa vara and a waddling gait. Radiography usually shows sclerosis of the ribs, flaring of the metaphyses, and a wide irregular growth plate, especially of the knees. A variant form of SMCD is spondylometaphyseal dysplasia Japanese type. It is characterized by spinal involvement comprising mild platyspondyly, vertebral body abnormalities, and end-plate irregularity., function:Type X collagen is a product of hyperthrophic chondrotocytes and has been localized to presumptive mineralization zones of hyaline cartilage., PTM:Prolines at the third position of the tripeptide repeating unit (G-X-Y) are

hydroxylated in some or all of the chains.,similarity:Contains 1 C

Background

This gene encodes the alpha chain of type X collagen, a short chain collagen expressed by hypertrophic chondrocytes during endochondral ossification. Unlike type VIII collagen, the other short chain collagen, type X collagen is a homotrimer. Mutations in this gene are associated with Schmid type metaphyseal chondrodysplasia (SMCD) and Japanese type spondylometaphyseal dysplasia (SMD). [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images