




GPX7 mouse mAb

Catalog No	YP-mAb-08317
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GPX7 GPX6 UNQ469/PRO828
Immunogen	Synthesized peptide derived from human GPX7 AA range: 77-127
Specificity	This antibody detects endogenous levels of GPX7 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse serum by affinity-chromatography using specific immunogen.
Dilution	WB 1 :500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Observed Band	
Calculated Molecular Weight	21kD
Cell Pathway	Secreted .
Tissue Specificity	Expressed in esophageal epithelial cells; expression is up-regulated after exposure to acidic bile acids.
Function	catalytic activity:2 glutathione + H(2)O(2) = glutathione disulfide + 2 H(2)O., sequence caution:Wrong choice of frame.,similarity:Belongs to the glutathione peroxidase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using GPX7 mouse mAb