



KIAA0182 Mouse mAb

Catalog No	YP-mAb-19516
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	GSE1 KIAA0182
Immunogen	Synthesized peptide derived from human KIAA0182. AA range:564-664
Specificity	This antibody detects endogenous levels of KIAA0182 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	GSE1 KIAA0182;Genetic suppressor element 1
Observed Band	
Calculated Molecular Weight	134kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

