



Glutaredoxin-1 Monoclonal Antibody

Catalog No	YP-mAb-10926
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GLRX GRX
Immunogen	Synthesized peptide derived from human Glutaredoxin-1 AA range: 8-58
Specificity	This antibody detects endogenous levels of human Glutaredoxin-1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, mouse,IgG
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Glutaredoxin 1;Glrx1;Glutaredoxin-1;Thioltransferase-1;TTase-1;GLRX GRX; Glutaredoxin-1 (Thioltransferase-1;TTase-1)
Observed Band	
Calculated Molecular Weight	12kDa
Cell Pathway	Cytoplasm.
Tissue Specificity	Blood,Brain,Lens,Liver,Skin,Spleen,Urinary bladder,
Function	function:Has a glutathione-disulfide oxidoreductase activity in the presence of NADPH and glutathione reductase. Reduces low molecular weight disulfides and proteins.,similarity:Belongs to the glutaredoxin family.,similarity:Contains 1 glutaredoxin domain.,
Background	This gene encodes a member of the glutaredoxin family. The encoded protein is a cytoplasmic enzyme catalyzing the reversible reduction of glutathione-protein mixed disulfides. This enzyme highly contributes to the antioxidant defense system. It is crucial for several signalling pathways by controlling the S-glutathionylation status of signalling mediators. It is involved in beta-amyloid toxicity and Alzheimer's disease. Multiple alternatively spliced transcript variants encoding the same protein have been identified. [provided by RefSeq,


Aug 2011],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using Glutaredoxin-1 Monoclonal Antibody