




TPD52 Mouse mAb

Catalog No	YP-mAb-00537
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Specificity	TPD52 Polyclonal Antibody detects endogenous levels of TPD52 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	Affinity purification
Dilution	WB: 1/500 - 1/2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	D52;hD52;N8L;PC 1;PrLZ;Protein N8;TPD52;tumor protein D52;PC-1
Observed Band	26kDa
Calculated Molecular Weight	24kDa
Cell Pathway	Mitochondrion .
Tissue Specificity	Mitochondrion .
Function	GTPase involved in the 5-carboxymethylaminomethyl modification (mnm(5)s(2) U34) of the wobble uridine base in mitochondrial tRNAs.
Background	Enables calcium ion binding activity and protein homodimerization activity. Involved in B cell differentiation. Located in endoplasmic reticulum and perinuclear region of cytoplasm.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of various lysates using TPD52 Mouse mAb (A10254) at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Mouse IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.