



LOXL1 Rabbit pAb

Catalog No	YP-Ab-18716
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	LOXL1 LOXL
Immunogen	Synthesized peptide derived from human LOXL1
Specificity	This antibody detects endogenous levels of LOXL1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LOL;LOXL;LOXL1;LOXL1 LOXL;Lysyl oxidase homolog 1 (Lysyl oxidase-like protein 1) (LOL)
Observed Band	80kDa/63kDa
Calculated Molecular Weight	63kDa
Cell Pathway	Secreted, extracellular space .
Tissue Specificity	Expressed in ocular tissues including the iris, ciliary body, lens and optic nerve. Not detected in the retina.
Function	Active on elastin and collagen substrates.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

