



MFRN1 mouse mAb

Catalog No	YP-mAb-10173
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLC25A37 MFRN MSCP HT015
Immunogen	Synthesized peptide derived from human MFRN1 AA range: 187-237
Specificity	This antibody detects endogenous levels of human MFRN1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MFRN;Mfrn1;Mitochondrial iron transporter 1;Mitoferrin 1;MSC;MSCP;SLC25A37;SLC25A37 MFRN MSCP HT015;Mitoferrin-1 (Mitochondrial iron transporter 1) (Mitochondrial solute carrier protein) (Solute carrier family 25 member 37)
Observed Band	37 kDa
Calculated Molecular Weight	
Cell Pathway	Mitochondrion inner membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

