



GATAD2A Mouse mAb

Catalog No	YP-mAb-19532
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	GATAD2A
Immunogen	Synthesized peptide derived from human GATAD2A. AA range:231-331
Specificity	This antibody detects endogenous levels of GATAD2A at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	GATAD2A;Transcriptional repressor p66-alpha (Hp66alpha) (GATA zinc finger domain-containing protein 2A)
Observed Band	
Calculated Molecular Weight	70kD
Cell Pathway	Nucleus speckle . Speckled nuclear localization requires both CR1 and CR2 regions. .
Tissue Specificity	Ubiquitous , both in fetal and adult tissues.
Function	Transcriptional repressor. Enhances MBD2-mediated repression. Efficient repression requires the presence of GATAD2B.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images