



Tubulin β mouse mAb(PT0354)

Catalog No	YP-Ab-04881
Isotype	IgG
Reactivity	Human
Applications	IHC-p;IF(paraffin section);WB
Gene Name	TUBB3 TUBB4
Protein Name	Tubulin β
Immunogen	Synthesized peptide derived from human Tubulin β
Specificity	This antibody detects endogenous levels of human Tubulin β
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:100-300, WB 1:500-2000
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	Tubulin beta-3 chain (Tubulin beta-4 chain;Tubulin beta-III)
Observed Band	42kD
Cell Pathway	Cytoplasm, cytoskeleton . Cell projection, growth cone . Cell projection, lamellipodium . Cell projection, filopodium .
Tissue Specificity	Expression is primarily restricted to central and peripheral nervous system. Greatly increased expression in most cancerous tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

