



PTMS Rabbit pAb

Catalog No	YP-Ab-19535
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	PTMS
Immunogen	Synthesized peptide derived from human PTMS. AA range:1-130
Specificity	This antibody detects endogenous levels of PTMS at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	PTMS;Parathymosin
Observed Band	
Calculated Molecular Weight	11kD
Cell Pathway	
Tissue Specificity	
Function	Parathymosin may mediate immune function by blocking the effect of prothymosin alpha which confers resistance to certain opportunistic infections.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

