



TMEM87A Mouse mAb

Catalog No	YP-mAb-19539
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	TMEM87A PSEC0094
Immunogen	Synthesized peptide derived from human TMEM87A. AA range:74-144
Specificity	This antibody detects endogenous levels of TMEM87A at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	TMEM87A PSEC0094;Transmembrane protein 87A
Observed Band	
Calculated Molecular Weight	61kD
Cell Pathway	Golgi apparatus membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	May be involved in retrograde transport from endosomes to the trans-Golgi network (TGN) .
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

