



COG1 Mouse mAb

Catalog No	YP-mAb-19540
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	COG1 KIAA1381 LDLB
Immunogen	Synthesized peptide derived from human COG1. AA range:661-731
Specificity	This antibody detects endogenous levels of COG1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	COG1 KIAA1381 LDLB;Conserved oligomeric Golgi complex subunit 1 (COG complex subunit 1) (Component of oligomeric Golgi complex 1)
Observed Band	
Calculated Molecular Weight	108kD
Cell Pathway	Golgi apparatus membrane ; Peripheral membrane protein ; Cytoplasmic side .
Tissue Specificity	
Function	Required for normal Golgi function.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

