




NUBP1 Mouse mAb

Catalog No	YP-mAb-19541
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive .
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	NBP; NBP1; CIAO5; NBP35; NUBP1
Observed Band	38kDa
Calculated Molecular Weight	35kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	NUBP1 is a member of the NUBP/MRP subfamily of ATP-binding proteins (Nakashima et al ., 1999 [PubMed 10486206]).
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HepG2 cells,
using NUBP1 Mouse mAb at 1:1000 dilution .
Secondary antibody: HRP-conjugated Goat anti-Mouse
IgG (H+L) at 1:10000 dilution .
Lysates/proteins: 25µg per lane .
Blocking buffer: 3% nonfat dry milk in TBST .
Detection: ECL Basic Kit (RM00020).
Exposure time: 1s .